

MONTHLY BULLETIN

ATTENTION: ALL FACILITIES

**UPDATES OCCURRING WITHIN THE CIS THE
WEEK OF - AUGUST 28TH, 2017**

The following solution areas have updates:


Lab - Pages 2 - 5

UPCOMING EVENTS


There are no upcoming events in CIS.

Lab Updates:

Changes in Heading Display Names: The heading for the Microbiology Referred Procedures is now clearly defined.



The Misc Blood Bank section heading now displays as Blood Bank Miscellaneous



Lab Updates:

Blood Bank Specimens on patients with known antibodies: Patients with known antibodies will now have a second set of labels (both Specimen label and Signature label) print out for the collection of a second tube to provide enough sample to complete the Antibody Investigation which is required for that patient. The Specimen Label ordered by the system will be labeled with ABID Tube, the Second Signature Label will remain the same. See example below. Like all Blood Bank Specimens, the labels need to be signed, dated and timed.

Exemptions are ordered placed on patients in the OR.

The diagram illustrates two sets of lab labels for a patient named Sally Jo. The top set, labeled 'Original labels', consists of a specimen label and a signature label. The specimen label includes the following information: 010200681, 09AUG17 1537, ZZ TEST, SALLY JO, DSCH, S M F, 17-221-00006A Blood, RT/RT, BB Order, 6.00mL Purple 6 BB, and Signature--:. The signature label includes: 010200681, 09AUG17 1537, ZZ TEST, SALLY JO, Collector's Signature:, and a certification statement: 17-221-00006A, I certify that I have checked the name and ID number of the patient and on the blood specimen label.

The bottom set, labeled 'additional labels ordered by the system for a patient with an antibody', is identical to the top set but includes 'ABID Tube' in the BB Order field, which is highlighted with a red box. The signature label is also identical to the top set.

Lab Updates:


Blood Bank Transfusion Alert for patients with hgb over 90g/L: This Alert will be activated when a crossmatch is being ordered on a patient whose hgb is over 90g/L (in the last 48 hrs). The Alert suggests that a transfusion is inappropriate based on the patient's most recent hgb, which is also contained in the Alert. Options allow canceling or continuing the crossmatch.

There is a Clinical Guidelines button at the bottom leading to a webpage from which this Alert is based.

If either of the following Red Cell Reasons are selected, the Alert will not fire:

1. Major Bleeding
2. Surgical Procedure likelihood of bleeding

The Alert will also not fire for orders placed through SurgiNet.



Lab Updates:

New Blood Bank Rule Alert: This Alert will be activated in PowerChart and FirstNet if a Group and Screen is being ordered on a patient and there still is a Group and Screen in date in the lab. No new specimen is required. The specimen in the lab is active and any red cells required can be matched on the specimen in the lab via a Telephone Request for Red Cells.

