Parent/ Child Orders


Parent/ Child orders exist when the orginal order has recurring orders attached.

For example:

The parent order of a timed study would be 'cbc w/auto diff daily-now for three times'. The children of this order would be the individual cbc orders for the three times

To view the child orders, a advanced fliter needs to be set:

1. Click on the elipsis button beside the display window



2. The Advance Filter window opens. Click on the box –Show individual instances for continuing order.

Advanced Filters		—
Display: Clinician Default		
Order s Active Window opens Vordered InProcess Future Incomplete Suspended	Inactive statuses: Inactive statuses: Inactive Inacti	Load the following types of orders Inpatient orders Ambulatory (In Office) orders Prescriptions Documented medications Miscellaneous
On Hold, Med Student Include all days Include number of days: 365	 Transfer/Canceled Include all days Include number of days: 	Show individual instances for continuing orders Inactive Plan Filter
Save	ave As Apply	New Remove Cancel

3. Click Apply.

On the orders page Parent /Child order will display with the black delta symbol *

		Ð		2016-Dec-21 11:32	Ordered	Blood Bank Comment	Blood, RT collect, Collect
		Click or	і віаск		(In-Lab)		
Distance		🔁 De	lta	 2017-Nov-08 09:28	Ordered	CBC w/ Auto Diff	Blood, Timed Study colle
Discharge							time(s)
	-	a		 			

When you click on the black delta the parent order will expand to display the child orders.

1	2017-Nov-08 09:28	Ordered	CBC w/ Auto Diff	Blood, Timed Study collect, Start: 2017-Nov-08 09:28, interval time(s)		
	2017-Nov-08 09:28	Ordere d (Sc	CBC w/ Auto Diff	Blood, Timed Study collect, Start: 2017-Nov-08 14:30		
	2017 Mar 08 09:28	Ordere dividual	CBC w/A Diff	Blood, Timed Study collect, Start: 2017-Nov-08 13:30		
	2017-N 09:28	Orde		Blood, Timed Study collect, Start: 2017-Nov-08 12:30		
	2017 - 107 - 00 09:28	d (Dis	CBC w/ AL Diff	Blood, Timed Study collect, Start: 2017-Nov-08 11:30		
	2017-Nov-08	Ordere	CBC w/ Auto Diff	Blood, Timed Study collect, Start: 2017-Nov-08 10:30		

The 'Show individual instances of continuing orders setting' enables the display of the 'parent' or initial order with a list of subsequent 'child' orders listed below for the next 12 hours..