



## How do I use my Genuair®?

1. Remove the cap (lightly squeeze the arrows).
2. Press the green button all the way down ONCE and release.
3. Check that the color control window is green. Green means ready.
4. Breathe out away from the mouthpiece.
5. Place your lips around the mouthpiece. Breathe in **strongly** and **deeply** through the mouthpiece. Keep breathing in even after you hear the inhaler “click”. Remove Genuair from mouth and hold breath for 5-10 seconds. Breathe out.
6. Check that the color control window has turned red.
7. If the window has not turned red repeat steps 4-6.
8. Replace cap.

*Note:* When a red band begins to appear in the dose counter this means you are nearing your last dose. The Genuair® locks after the last dose.

## Care of a Genuair®

1. Wipe the mouthpiece with a dry tissue or cloth.



For video instruction:  
[lungsask.ca/inhaler-resources](https://lungsask.ca/inhaler-resources)