

## How do I use my Nebulizer<sup>®</sup>?

1. Ensure nebulizer is plugged in & functioning.
2. Attach tubing to compressor air outlet.
3. Unscrew the top of the medicine cup.
4. Add the medication dose to the cup.
5. Screw on the top.
6. Turn on the compressor; ensure medication is misting.
7. Place either a mask on face or mouthpiece in mouth.
8. Complete the treatment.
9. Turn off the compressor.

## Care of a Nebulizer<sup>®</sup>

1. Wash mask and nebulizer medication chamber in warm, soapy water.
2. Rinse well and allow them to air-dry before reuse.

For video instruction:  
[lungsask.ca/inhaler-resources](http://lungsask.ca/inhaler-resources)